



**Model No.: 74-4122
Tonsillectomy Set**


1



Yankaur Suction Tube, SS

- Qty: 1pc

2



Davis Boyle Mouth Gag, SS

- 5 Blades
- Qty: 1pc

3



Knife Handle, SS

- Size: No. 7
- Qty: 1pc

4



Scissors, Metzenbaum, TC

- Curved
- Size: 17.8cm
- Qty: 1pc

5



Forceps, Mikulicz Tonsil, SS

- Size: 19cm
- Qty: 1pc

6



Towel Clamp, Backhaus, SS

- Size: 13cm
- Qty: 4pcs

7



Forceps, Ballenger Sponge, SS

- Straight, Serrated
- Size: 17.8cm
- Qty: 1pc

8&9



Needle Holder, Crile Wood, TC

- Size: 15cm & 20cm
- Qty: 1pc each

10,11&12




Barnhill Adenoid Curette, SS

- No. 0, 1, 2
- Size: 21.6cm
- Qty: 1pc each



**Model No.: 74-4122
Tonsillectomy Set**

13



Retractor, Love Uvula, SS

- Size: 18mm
- Qty: 1pc

14



Meltzer Adenoid Punch, SS

- Triangular Jaws
- Size: 1
- Qty: 1pc

15



Wieder Tongue Depressor SS

- Size: Large
- Qty: 1pc


16



Forceps, White Tonsil, SS

- Curved
- Size: 19cm
- Qty: 1pc


17



Forceps, Schindt Tonsil, SS

- Size: 19cm
- Qty: 6pcs


18



Forceps, Sawtell Tonsil, SS

- Size: 19cm
- Qty: 1pc


19



**Fisher Tonsil Dissector & Pillar
Retractor, SS**

- Size: 21.6cm
- Qty: 1pc

20



**Hurd Tonsil Knife & Dissector and
Pillar Retractor, SS**

- Size: 22.9cm
- Qty: 1pc

21



Tonsil Snare, SS

- Qty: 1pc