

**Model No.: 74-4091
Cesarean Section Set**

1



Retractors, Richardson, SS

- Grip Handle
- Length: 26.5cm,
- Blade Size: 3.8 cm x 2.5 cm
- Qty: 1pc

2



Retractors, Richardson, SS

- Grip Handle
- Length: 26.5cm,
- Blade Size: 3.8 cm x 2.5 cm
- Qty: 1pc

3



Retractors, Deaver, SS

- Blade Size: 7.5cm Wide
- Size: 30cm
- Qty: 1pc


4



Retractor, De Lee Universal, SS

- Rounded Jaws
- Blade: 7cm
- Wide: 5.1cm
- Size: 24 cm
- Qty: 1pc


5



Forceps, Tissue, SS

- Teeth: 4x5
- Size: 14cm
- Qty: 2pcs

6



Forceps, Adson, SS

- Teeth: 1x2 with Serrations
- Delicate
- Size: 12cm
- Qty: 2pcs

7



Forceps, Tissue, Russian, SS

- Size: 20cm
- Qty: 1pc

8



Scissors, Dissecting, Vital Mayo, TC

- Straight
- Size: 17cm
- Qty: 2pcs

9



Scissors, Dissecting, Vital Mayo, TC

- Curved
- Size: 17cm
- Qty: 2pcs



Model No.: 74-4091
Cesarean Section Set

10



Scissors, Dissecting, Vital Metzenbaum, TC

- Delicate
- Blunt/Blunt, Straight
- Size: 18cm
- Qty: 2pcs

11



Scissors, Umbilical, SS

- Straight
- American Pattern
- Size: 10cm
- Qty: 1pc


12



Scissors, Lister Bandage, SS

- Angled, Probe Point
- Size: 14cm
- Qty: 2pcs

13



Forceps, Towel Roeder, SS

- Ball Stops
- Size: 14cm
- Qty: 4pcs


14



Forceps, Crile, SS

- Regular Pattern, Curved
- Size: 14cm
- Qty: 12pcs

15



Forceps, Crile, SS

- Regular Pattern, Straight
- Size: 14cm
- Qty: 6pcs


16



Forceps, Tissue Holding, Babcock, SS

- Size: 16cm
- Qty: 2pcs

17



Forceps Tissue, Allis, SS

- Teeth: 4x5
- Size: 15cm
- Qty: 6pcs

18



Forceps, Artery Pean, SS

- Straight
- Size: 16.5cm
- Qty: 6pcs



Model No.: 74-4091
Cesarean Section Set


19



Forceps Sponge Holding, Ballenger, SS

- Straight, Serrated Loop Jaws
- Size: 18cm
- Qty: 4pcs


20



Forceps Tissue Holding, Pennington, SS

- Triangular serrated jaws
- Size: 20cm
- Qty: 4pcs

21



Needle Holders, Mayo Hegar, SS

- Straight, Titanium Delicate Touch
- Dust Jaw Surface Foruse With 3-0 Suture & Smaller
- Size: 16cm
- Qty: 2pcs

22



Suction Tubes, Abdominal, Poole, SS

- Straight, Mirror Finish, 30FR
- Size: 24cm
- Qty: 1pc

23



Suction Tubes, Yankauer, SS

- Half Curved
- With additional slight curved near the tip, Large, light, hollow handle.
- Mirror Finish
- Size: 27cm
- Qty: 1pc

24



Handle Knife, SS

- Handle No. 3
- Mirror Finish
- Blade No.: 10, 11, 12 and 15
- Size: 12.5cm
- Qty: 2pcs

25



Handle Knife, SS

- Handle No. 4
- Mirror Finish
- Blade No.: 20, 21, 22, 23, 24 and 25
- Size: 13.5cm
- Qty: 1pc

26&28



Forceps Artery, Rochester, Ochsner (Kocher), SS

- Teeth: 1x2, Serrated Jaws
- Straight
- Size: 20cm
- Qty: 2pcs each

27&29



Forceps Artery, Rochester, Ochsner (Kocher), SS

- Teeth: 1x2, Serrated Jaws
- Curved
- Size: 20cm
- Qty: 2pcs each



Model No.: 74-4091
Cesarean Section Set

30



**Forceps, Heaney Ballentine,
Hemostatic, SS**

- Single Tooth
- Curved Heavy
- Size: 20cm
- Qty: 8pcs

35



Instrument Tray, SS

- Qty: 1pc

31&32



Forceps, Artery, Kelly, SS

- Straight
- Size: 14cm & 16cm
- Qty: 6pcs each

33&34



Forceps, Artery, Kelly's, SS

- Curved
- Size: 14cm & 16cm
- Qty: 6pcs each

The Complete Kit Packed in a Carry Bag