




Model No.: 74-4051
Haemorrhoidectomy Set


1



Forceps, Forester Sponge Holding, SS

- Size: 18cm
- Serrated
- Qty: 2pcs

2



Forceps, Towel, Backhaus, SS

- Size: 9cm
- Qty: 4pcs

3



Scalpel Handle, SS

- Size; No. 3
- Qty: 1pc

4



Scissors, Mayo Dissection, SS

- Curved
- Size: 14.5cm
- Qty: 1pc

5



Forceps, Tissue Gillies, SS

- Size: 15cm
- Teeth: 1x2
- Qty: 1pc

6



Forceps, Mc Givney, SS

- Size: 19 cm
- Qty: 1pc

7



Forceps, Leriche Haemostatic, SS

- Straight
- Teeth: 1x2
- Size: 15cm
- Qty: 4pcs

8



Sims Rectal Specula, SS

- Blade: 9cm
- Size: 15cm
- Qty: 1pc

9




Kelly Sphincteroscope With Obturator, SS

- Diameter: 27x50mm
- Qty: 1pc



Model No.: 74-4051
Haemorrhoidectomy Set

10



Brodie Probe, SS

- Size: 20cm
- Qty: 1pc


11



Fergusson Angiotribe, SS

- Curved
- Size: 20cm
- Qty: 1pc

12



MC Givney, Haemorrhoidal Ligator Complete, SS

- Qty: 1pc

13



Ligator Rings, SS

- Pack of 100
- Qty: 1pc

14



Needle Holder, Mayo Hegar, TC

- Size: 16cm
- Qty: 1pc

15



Gallipot, SS

- Diameter: 10cm
- Qty: 2pcs

16



Kidney Dish, SS

- Size: 25cm
- Qty: 1pc

17



Sterilization Container, SS

- Size: 28.5x28x10 cm

18



Wire Mesh Basket, SS

- Qty: 1pc



Model No.: 74-4051
Haemorrhoidectomy Set

19



Identification Labels

- Color: red
- Qty: 2pcs

The Complete Kit Packed in a Carry Bag