



Model No.: 74-4087
Head & Neck Set

1



Scalpel Handle, SS

- Size: No. 3
- Qty: 1pc

2



Scalpel Handle, SS

- With Graduation cm/inch No. 3
- Size: No. 1
- Qty: 1pc


3



Forceps, Halsted Mosquito, SS

- Curved
- Qty: 2pcs

4



Forceps, Halsted Mosquito, SS

- Straight
- Qty: 16pcs

5



Forceps, Jacobson Mosquito, SS

- Qty: 2pcs

6



Forceps, Jacobson Hemostatic, SS

- Qty: 4pcs


7



Forceps, Artery Kelly, SS

- Qty: 4pcs


8



Forceps, Artery Crile, SS

- Qty: 6pcs

9



Forceps, Tissue Allis, SS

- Qty: 2pcs



Model No.: 74-4087
Head & Neck Set

10



Forceps, Tissue Holding, Babcock, SS

- Qty: 2pcs


11



Forceps, Artery Pean, SS

- Qty: 2pcs


12



Mixer Micro Line, SS

- Size: Pediatric
- Qty: 1pc

13



Forceps, Towel Backhaus, SS

- Qty: 4pcs

14



Forceps, Lahey Guiter Vulsellum, SS

- Qty: 4pcs


15



Mccabe Facial Nerve Dissector, SS

- Qty: 1pc

16



Needle Holder, Halsey Eye, SS

- Qty: 2pcs

17



Needle Holder, Derf Eye, SS

- Qty: 1pc

18



Needle Holder, Mayo Hegar, SS

- Qty: 1pc



Model No.: 74-4087
Head & Neck Set


19



Scissors, Iris, SS

- Qty: 1pc

20



Scissors, Knapp Strabismus, SS

- Qty: 1pc

21



Scissors, Metzenbaum, Dissecting, SS

- Curved
- Qty: 1pc

22



Scissors, Metzenbaum Lipsett, SS

- Qty: 1pc

23



Scissors, Metzenbaum, Dissecting, SS

- Straight
- Qty: 1pc

24



Scissors, Kaye Face Lift, SS

- Qty: 1pc

25



Scissors, Potts Tenotomy, SS

- Qty: 1pc

26



Scissors, Mayo, Dissecting, SS

- Straight
- Qty: 1pc

27



Scissors, Mayo Dissecting, SS

- Curved
- Qty: 1pc



Model No.: 74-4087
Head & Neck Set

28



Retractor, Weitlaner, SS

- Qty: 1pc


29



Retractor, Baby Weitlaner, SS

- Qty: 1pc

30



Retractor, Perineal Gelpi, SS

- Qty: 2pcs


31



Retractor, USA, SS

- Qty: 1pc

32



Retractor, Green Goiter, SS

- Qty: 2pcs

33



Weder Tongue Depressor, SS

- Qty: 1pc

34



Retractor, Richardson, SS

- Qty: 2pcs


35



Retractor, Sistrunk, SS

- Qty: 2pcs

36



Retractor, Senn, SS

- Qty: 2pcs



Model No.: 74-4087
Head & Neck Set


37



Forceps, Adson, SS

- Qty: 2pcs

38



Forceps, Adson Brown Tissue, SS

- Qty: 2pcs

39



Forceps, Cushing Tissue, SS

- Qty: 2pcs


40



Forceps, Debaquey Diethrich Tissue, SS

- Qty: 2pcs


41



Forceps, DeBauey Tissue, SS

- Qty: 2pcs


42



Forceps, Martin Tissue, SS

- Qty: 2pcs

43



Trousseau Jackson Tracheal, SS

- Qty: 1pc


44



Ferris Smith Sewall Orbital, SS

- Different Size
- Qty: 1pc each

45




Suction Tube, Frazier, SS

- Qty: 1pc



Model No.: 74-4087
Head & Neck Set


46



Tube Suction, Frazier, SS

- Qty: 2pcs


47



Double Prong Micro Hook, SS

- Qty: 3pcs


48



Jackson Tracheal Tenaculum, SS

- Qty: 2pcs

49



Nasal Tenacula, SS

- Qty: 2pcs


50



Freer Septum Elevator, SS

- Qty: 1pc

51



Krayenbuehi Hook, SS

- Qty: 1pc

52



Instrument Tray, SS

- Qty: 1pc

The Complete Kit Packed in a Carry Bag