



Model No.: 74-4098
Lacrimal Set

1



Clamp Towel, Backhaus, SS

- Qty: 1pc

2



Bishop Harman A/C Irrigator

- Cannula 19G, Silicone Bulb & Adapter
- Qty: 1pc


3



Beyer, Rongeur, SS

- Single Action
- Qty: 1pc

4



Retractor Stevenson, SS

- 3 x 3 Prongs
- Qty: 1pc

5



Kerrison Bone Nibbling Rongeur, SS

- Wide: 4.5mm
- Qty: 1pc


6



Nasal Speculum, SS

- Size: Adult
- Qty: 1pc

7



Castroviejo Lacrimal Dialator, SS

- Qty: 1pc

8



Mallet

- Qty: 1pc

9



Free Periosteal Elevator, SS

- Double End
- Qty: 1pc



Model No.: 74-4098
Lacrimal Set

10



Forceps, Sauer Fixation, SS

- Size: 1x2
- Qty: 1pc

11



Forceps Hartman Mosquito, SS

- Straight
- Size: 10cm
- Qty: 2pcs

12



Forceps Hartman Mosquito, SS

- Curved
- Size: 10cm
- Qty: 2pcs


13



Crawford Retrieving Hook, SS

- Qty: 1pc

14



Lacrimal Cannula, SS

- Malleable Tip
- Straight
- Size: 23G
- Qty: 1pc

15



Lacrimal Cannula, SS

- Curved Gently
- Size: 23G
- Qty: 1pc


16



West Lacrimal Chisel, SS

- Qty: 1pc


17



Pigtail Probe, SS

- Qty: 1pc

18




Freer Lacrimal Chisel, SS

- Curved
- Qty: 1pc



Model No.: 74-4098
Lacrimal Set


19



Crawford Lacrimal Intubation Set, SS

- 2 probes with silicon tubing
- Qty: 1pc

20



Bowman Lacrimal Probe Set, SS

- Qty: 1pc

21



Forceps, Dressing, SS

- Serrated Delicate
- Straight
- Size: 14cm
- Qty: 1pc


22



Forceps Bonaccolto Conjunctiva, SS

- Qty: 1pc

23



Forceps, Castroviejo Suturing, SS

- Straight with tying format
- Teeth: 1x2
- Tip Size: 0.12mm
- Qty: 1pc

24



Retractor, Knapp, SS

- 4 Prong, Blunt
- Qty: 2pcs

25



Scissors, Stevens Tenotomy, SS

- Curved
- Qty: 2pcs

26



Scissors, Eye, SS

- Straight, Standard
- Qty: 2pcs

27



Needle Holder, Castroviejo, SS

- Delicate Jaws
- Curved without Lock
- Qty: 1pc



Model No.: 74-4098
Lacral Set

28



Scissors, McPherson Westcott Stitch, SS

- Sharp Tip
- Qty: 2pcs

29



Bard Parker Blade Handle, SS

- Double
- Qty: 1pc

30



Silicone Mat

Sterilization Box, SS

- With Two Silicon Mats
- Size: Large
- Qty: 1pc

The Complete Kit Packed in a Carry Bag