



Model No.: 74-4110
Nephrotomy Set

1



Instrument Tray, SS

- With Cover
- Qty: 1pc

2



Wire Mesh, SS

- Size: 48x24x5cm
- Qty: 1pc

3



Scalpel Handle, SS

- No.: 3L
- Qty: 1pc

4



Scissors, Metzenbaum, SS

- Curved
- Size: 18cm
- Qty: 1pc

5



Forceps, Dressing, SS

- Standard, Straight
- Size: 14.5cm
- Qty: 2pcs

6



Forceps, Tissue, SS

- Standard, Straight
- Teeth: 1 x 2
- Size: 20cm
- Qty: 1pc

7



Forceps, Tissue, Allis, SS

- Size: 15cm
- Qty: 2pcs

8



Forceps, Artery, Halstaed Mosquito, SS

- Straight
- Size: 12.5cm
- Qty: 2pcs

9



Forceps, Artery, Halstaed Mosquito, SS

- Curved
- Size: 12.5cm
- Qty: 10pcs



Model No.: 74-4110
Nephrotomy Set

10



Forceps, Artery, Kocher, SS

- Straight
- Teeth: 1 x 2
- Size: 16cm
- Qty: 10pcs

11



Retractor, Richardson, SS

- Blade Size: 28 x 20 mm
- Size: 24cm
- Qty: 4pcs

12



Retractor, Richardson, SS

- Blade Size: 36 x 28mm
- Size: 24cm
- Qty: 2pcs

13



Retractor Set, Roux, SS

- Qty: 2pcs

14



Retractor, Deaver, SS

- Size: 30cm, 50mm(Width)
- Qty: 1pc

15



Retractor, Volkmann, SS

- Semi-Sharp, 4 prongs
- Size: 21.5cm
- Qty: 2pcs


16



Retractor, Volkmann, SS

- Semi-Sharp, 6 Prongs
- Size: 21.5cm
- Qty: 2pcs

17



Retractor, Cushing, SS

- Size: 20cm
- Width: 10mm
- Qty: 10pcs

18



Retractor, US Army, SS

- Set of 2
- Qty: 1pc



Model No.: 74-4110
Nephrotomy Set


19



Spatula, Abdominal, SS

- Malleable
- Size: 33cm
- Wide: 30mm
- Qty: 1pc

20



Probe with Eye, SS

- Size: 13cm x 2mm
- Qty: 1pc


21



Conductor, Ligature, König, SS

- Size: 19.5cm x 3mm
- Qty: 1pc


22



Needle, ligature, Deschamps, SS

- Left Blunt
- Size: Medium, 20cm
- Qty: 1pc

23



Needle, ligature, Deschamps, SS

- Right Blunt
- Size: Medium, 20cm
- Qty: 1pc

24



Needle Holder, Mayo Hegar, SS

- Size: 16cm
- Qty: 1pc


25



Needle Holder, Mayo Hegar, SS

- Size: 18cm
- Qty: 2pcs

26



Scissors, Nelson (Metzenbaum), SS

- Curved
- Size: 23cm
- Qty: 2pcs

27



Forceps, Dressing, Potts Smith, SS

- Straight
- Size: 21cm
- Qty: 2pcs



Model No.: 74-4110
Nephrotomy Set

28



Forceps, Dressing, Potts Smith, SS

- Size: 25cm
- Straight
- Qty: 1pc

29



Needle Holder, Mayo Hegar, SS

- Size: 24cm
- Qty: 1pc


30



Forceps, Artery, Kocher, SS

- Teeth: 1x2
- Straight
- Size: 16cm
- Qty: 1pc


31



Forceps, Artery, Mixter, SS

- Curved
- Size: 19cm
- Qty: 8pcs


32



Forceps, Intestinal, Tissue, Allis, SS

- Size: 25cm
- Qty: 2pcs


33



Forceps, Kidney Pedicle, Guyon, SS

- Size: 24cm
- Qty: 4pcs

34



Forceps, Kidney Pedicle, Herrick, SS

- Size: 23cm
- Qty: 2pcs


35



Forceps, Clamp, Wertheim Cullen, SS

- Size: 21.5cm
- Qty: 2pcs

36



Forceps, Gall Duct, Lahey, SS

- Size: 23cm
- Qty: 2pcs



Model No.: 74-4110
Nephrotomy Set

37



Retractor Richardson, SS

- Blade Size: 52 x 22 mm
- Size: 24cm
- Qty: 2pcs

38



Retractor Richardson, SS

- Blade Size: 65 x 50 mm
- Size: 26cm
- Qty: 2pcs

39



Needle Holder, Masson, SS

- Size: 27cm
- Qty: 2pcs

40



Rongeur, Bone, Stille Luer, SS

- Curved
- Size: 22cm
- Qty: 1pc

41



Raspatory, Rib, Doyen, SS

- Left
- Size: Adult, 17cm
- Qty: 1pc

42



Raspatory, Rib, Doyen, SS

- Right
- Size: Adult, 17cm
- Qty: 1pc


43



Raspatory, Alexander, SS

- Size: 20cm
- Qty: 1pc

44



Forceps, Bone Cutting, Horsley, SS

- Size: 27cm
- Qty: 1pc

45



Forceps, Bone Holding, Semb, SS

- With Ratchet
- Size: 20cm
- Qty: 1pc



Model No.: 74-4110
Nephrotomy Set

46



Shears, Rib, Giertz Stille, SS

- Size: 27cm
- Qty: 1pc

47



Forceps, Dressing, SS

- Standard
- Size: 25cm
- Qty: 1pc

48



Forceps, Tissue, SS

- Standard
- Teeth: 1x2
- Size: 25cm
- Qty: 2pcs


49



Forceps, Tissue, Potts Smith, SS

- Straight
- Size: 25cm
- Qty: 2pcs


50



Scissors, Mayo, SS

- Straight
- Size: 23cm
- Qty: 2pcs


51



Scissors, Mayo, SS

- Curved
- Size: 23cm
- Qty: 1pc

52



Pin, Instrument Holder, Mayo, SS

- Size: 14cm
- Qty: 2pcs

The Complete Kit Packed in a Carry Bag