



Model No.: 74-4120
Transphenoidal Set

1



Stille Chisel, SS

- Size: 20.5cm (Length), 10mm Tip
- Qty: 1pc

2



Forceps, Gross Maier, SS

- Straight
- Size: 20.5cm
- Qty: 1pc

3



Hooklet Cushing, SS

- Wide: 10mm
- Size: 24cm
- Qty: 1pc

4



Forceps Tampon Ulrich, SS

- Straight
- Size: 22.5cm
- Qty: 1pc

5



Scissors, Strabismus, SS

- Curved
- Blunt/Blunt
- Size: 11.5 cm
- Qty: 1pc

6



Forceps, Dressing, Potts-Smith, TC

- Size: 18.5cm
- Qty: 1pc

7



Forceps, Ear Dressing, Jansen, Bayonet, TC

- Size: 21.5cm
- Qty: 1pc


8



Forceps, Gillies Dissecting, SS

- 1x2mm Tip
- Size: 15cm
- Qty: 1pc

9



Forceps, Jansen Middleton Nasal Cutting, SS

- Size: 19cm
- Qty: 1pc



Model No.: 74-4120
Transphenoidal Set

10



Nasal Speculum, SS

- Cottle No. 4
- Size: 85mm (Length)
- Qty: 1pc

11



**Pituitary Specula, Papavero
Caspar, SS**

- Size: 80x13mm
- Qty: 1pc

12



**Pituitary Specula, Papavero
Caspar, SS**

- Size: 90x13mm
- Qty: 1pc

13



Kerrison Punch 40° Forward, SS

- 1mm Bite
- Size: 20 cm (Length)
- Qty: 1pc

14



Kerrison Punch 40° Forward, SS

- 2mm Bite
- Size: 20 cm (Length)
- Qty: 1pc


15



Cushing Pituitary Rongeurs, TC

- 2x10mm Tip
- Straight
- Size: 20cm
- Qty: 1pc


16



Cushing Pituitary Rongeurs, TC

- 2x10mm Tip
- Up
- Size: 20cm
- Qty: 1pc

17



Cushing Pituitary Rongeurs, TC

- 2x10mm Tip
- Down
- Size: 20cm
- Qty: 1pc

18




Forceps, Lanholt Hypophyseal, SS

- Size: 20cm
- Qty: 1pc



Model No.: 74-4120 Transphenoidal Set


19



Pituitary Ring Curette, SS

- 4mm Cutting Ring
- Straight
- Size: 23cm
- Qty: 1pc


20



Pituitary Ring Curette, SS

- 4mm Cutting Ring
- Cutting To Right
- Size: 23cm
- Qty: 1pc


21



Pituitary Ring Curette, SS

- 4mm Cutting Ring
- Cutting To Left
- Size: 23cm
- Qty: 1pc


22



Pituitary Ring Curette, SS

- 4mm Cutting Ring
- Angled Down
- Size: 23cm
- Qty: 1pc

23



Pituitary Ring Curette, SS

- 4mm Cutting Ring
- Angled Up
- Size: 23cm
- Qty: 1pc


24



Sickle Knife, Mikael, SS

- Size: 26cm
- Qty: 1pc


25,26,27,28&29



Cushing Pituitary Spoon, SS

- Length: 23.5cm
- Wide: 5mm, 6mm, 8mm, 10mm & 12mm
- Size: 1, 2, 3, 4 & 5
- Qty: 1pc each


30



Insulated Yasargil Bayonetted Bipolar, SS

- 1.0mm Tip Tapered
- Reusable cable
- Size: 26cm
- Qty: 1pc

31



Fergusson Neurosuction Tubes, SS

- Size: 2.5mm (W) x 110mm (L)
- Qty: 1pc



Model No.: 74-4120
Transphenoidal Set


32



Fergusson Neurosuction Tubes, SS

- Size: 3.0mm (W) x 110mm (L)
- Qty: 1pc

33



Fergusson Neurosuction Tubes, SS

- Size: 4.0mm (W) x 110mm (L)
- Qty: 1pc

34



Instrument Tray, SS

- Qty: 1pc

The Complete Kit Packed in a Carry Bag