



Model No.: 74-4126
Varicose Vein Stripping Set

1



Instruments Box, SS

- Size: 32x15x6cm
- Qty: 1pc

2



Scissors, Iris, SS

- Curved, Blunt-Blunt
- 25mm Blades
- Size: 11.4cm
- Qty: 1pc

3



Scissors, Mayo Stille Operating, SS

- Straight
- Size: 17cm
- Qty: 1pc

4



Scissors, Metzenbaum-Nelson Dissecting, SS

- Blunt/Blunt, Curved
- Size: 18cm
- Qty: 1pc

5



Scalpel Handle, SS

- Standard, Solid
- Size: No. 4
- Qty: 1pc

6



Forceps, Micro Mosquito Artery, SS

- Straight
- Size: 12.5cm
- Qty: 4pcs

7



Forceps, Pean Artery, SS

- Straight
- Delicate Pattern
- Size: 14cm
- Qty: 2pc

8



Forceps, Pean Artery, SS

- Curved
- Delicate Pattern
- Size: 14cm
- Qty: 2pc

9




Forceps, Kocher, Ochsner, Artery, SS

- 1x2 teeth, Curved
- Size: 16cm
- Qty: 2pcs



Model No.: 74-4126 Varicose Vein Stripping Set

10



Nabatoff Vein Probe Complete Set, SS

- Size: 138x100x25mm
- Qty: 1pc

11



Forceps, Standard Dressing, SS

- Narrow
- Size: 14cm
- Qty: 1pc


12



Forceps, Tissue, SS

- Teeth: 1x2
- Size: 14cm
- Qty: 1pc

13



Forceps, Semken Dressing, SS

- Size: 12.5cm
- Qty: 1pc

14



Forceps, Dressing, SS

- With Dissector End
- Size: 15.5cm
- Qty: 1pc

15



Retractor, Roux, SS

- Set of 3
- Qty: 1pc

16



Mayo-Hegar Needle Holder, TC

- Size: 18cm
- Qty: 1pc

17



Ryder Needle Holders, TC

- Size: 15cm
- Qty: 1pc

18



Forceps, Pean Murphy, SS

- Size: 14cm
- Qty: 1pc

The Complete Kit Packed in a Carry Bag